RADIATION COUNTER TUBE, self quenching, with mica window, for the detection of alpha and beta particles and photons down to very low energies.

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APPLICATION: amongst others for laboratory, medical and industrial tracer technics, for thickness measuring instruments, and for laboratory or general survey instruments.

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
min. 10^{10} counts
comptes
Zählungen

EFFICIENCY. If exposed to a radium source of $2.39 \pm 1\%$ millicurie the counting rate is 7210 counts/min for 2.5 m röntgen/hour.

For beta and alpha particles the counting rate will, of course, depend on the type of source and its geometry and self-absorption. Beta particles arriving at the window with an energy of more than 50 keV will be counted.

For alpha particles the window stopping power is equivalent to 3 cm of air and particles arriving at the window with an energy of approximately 5 MeV will be counted.

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A rubber cap, easily detached or slipped over the end, is provided with each counter to protect the mica window if the counter is not used for the measurement of low energy radiation.

The counter is delivered in an aluminium container, connected to the counter cathode and which need not be removed for testing the tube. It is possible to use the counter while still inside the container for measuring gamma radiation and high energy particles.

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